
6/95-Pr/Sf-

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$$f = \frac{S \cdot d \cdot b^2}{2 m_1 \cdot t^2 \cdot L} \quad (II)$$

The device can be used to demonstrate the gravitational attractive force between masses and to determine the gravitational constant  $f$ .

#### Measurement principle (Fig. 1)

The static equilibrium (position I) of a torsion pendulum with a period of approx. 10 min. is disturbed by a change in position of the outer masses  $m_1$ , which affect the dumbbell-shaped pendulum body. The oscillations become damped and the pendulum takes up a new equilibrium position (position 2). The angle between the two positions is a measure of the active gravitational force.

The oscillation of the pendulum, which is equipped with a concave mirror, is indicated by a light pointer. This is possible using

- visible light directly on a mm-scale or
- infrared light along with the IR position detector (332 11). The latter enables the measurement values to be plotted or evaluated by computer.

The gravitational constant  $f$  can be obtained from the oscillation curve with respect to time, the mass  $m_1$  and the geometry of the arrangement using either the end deflection method or (in a quicker process) the acceleration method.

In the *end deflection method* the torsion pendulum period  $T$  and the distance  $S$  separating the light pointer positions are evaluated for the two equilibrium positions:

$$f = \frac{\pi^2 \cdot b^2 \cdot d \cdot S}{m_1 \cdot T^2 \cdot L} \quad (I)$$

In the *acceleration method*, the acceleration of the torsion pendulum  $a = 2s/t^2$  is evaluated using the masses  $m_1$  after the equilibrium position has been disturbed:

$$f = \frac{S \cdot d \cdot b^2}{2 m_1 \cdot t^2 \cdot L} \quad (II)$$