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3.2 Mounting the stand base on the wall (required for wall mounting as shown in Fig. 3 before using the setup for the first time.)

Important:

Wall mounting of the stand material is not suitable for the experiment setup with the IR position detector 332 11 (see section 3.1.2; Fig. 4).

Correctly mount the stand base on its apex hole as shown in Fig. 5.1 - 3.3 and align parallel to the wall using levelling screws; use solid base material if the levelling screws press into the wall.

3.3 Preparing for the experiment

Assemble the arrangement - initially without the large lead balls - on the wall as shown in Fig. 3 (or similarly on the optical bench on a sturdy lab bench).

Align the lamp housing so that the reflected beam is projected unobstructed onto the scale, which should be positioned at least 5 m distant.

Align the filament so that it is vertical before attaching the condenser with screen holder and slit (vertically!) to the lamp housing.

Move the insert in the lamp housing so that a focused image of the filament is projected onto the torsion balance mirror. To do this, hold a white sheet of paper, on which the filament forms an image, directly in front of the balance.

Move the lamp housing so that the slit is focused onto the scale.

If necessary, adjust the zero-point as described in section 3.5.

Attach the lead balls and bring the arrangement to an extreme position. In carrying out this step make absolutely sure that the housing does not come into contact with either your finger or a lead ball.

3.4 Experiment procedure; measuring example

Important!

After the experiment has been set up as described in section 3.3, it must be left undisturbed for at least two hours so that the pendulum can settle in its equilibrium state.

Turn the ball carrier ⑨ without touching the housing.

To shorten this settling time, make use of the diamagnetism of the lead: each time one of the pendulum lead balls approaches the glass cover, hold up a strong magnet - without touching the glass - until the motion of the system is reversed.

The experiment setup must not be subjected to any shocks or temperature variations which may cause convention in the torsion balance housing while recording measured values (see Fig. 6)

Check the zero-point stability before starting the measurement. If necessary, observe and note the zero-point fluctuations for at least 10 minutes so that you can determine an average value for x_0 .

At time $t = 0$ rotate the carrier with the lead balls quickly from one extreme position to the other. However, in carrying out this movement, make sure that fingers or lead balls do not come into contact with the housing. Start the stop-clock immediately after shifting the carrier.

Read off and note the position of the light pointer on the scale every 30 s (Fig. 6) for at least 3 periods when using the end displacement method and for 1 period when using the acceleration method.