Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

10/1989

Gebrauchsanweisung
Instruction Sheet
Instrucciones de Servicio

558 79

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Determining Planck's constant by means of the photoelectric effect using the compact arrangement

In the case of the photoelectric effect, the energy of the released electrons depends only on the frequency ν of the incident light. The quantum of action h is determined using this functional relationship.

Electrons may be knocked out from a metal surface due to irradiation with light (photoelectric effect). The number of electrons depends on the intensity of the incident light, but their energy, however, only on its frequency ν . This situation may be clarified if light is considered as a stream of "photons" and it is assumed that each photoelectron is released by an individual photon possessing the energy $h\nu$.

If we designate the work function of the electron out of the photocathode as A_C , we obtain the "Einstein relationship" for its kinetic energy.

$$\frac{1}{2} m v_0^2 = h\nu - A_C \quad (1)$$

m_0 : rest mass of the electron

A_C : work function dependent on the cathode material.

The measurements will be carried out using an alkaline photocell:

Light strikes a layer of potassium after passing through a ring-shaped anode (platinum wire). The photoelectrons reach the anode and are measured as photocurrent I_p . The photocurrent decreases if the electrons are allowed to come up against a negative anode voltage U (Fig. 1).


Fig. 1: Setup for the photoelectric effect, schematic.

The kinetic energy of the electrons varies, so that less and less electrons can overcome the retarding potential for increasing U , dependent on the state of the conduction band from which the electrons originate, the magnitude of their velocity component in the direction of the electric retarding field and on the collision processes they carry out on their path. $I_p = 0$

when even the electrons from the upper edge of the conduction band with emission velocities in the direction of the electric retarding field can no longer reach the anode.

The required negative anode potential to achieve this state is designated by $U_0(\nu)$.

In accordance with (1):

$$h \cdot \nu - A = e \cdot U_0 \quad (2)$$

A : work function

(see note 3 at the end of the section)

If the frequency of the incident light is increased by $\Delta\nu$, the electron energy increases by $h\Delta\nu$, so that the negative anode potential must be increased ΔU_0 in order to bring the current back to zero. Here, the following equation applies

$$e\Delta U_0 = h\Delta\nu \quad (e = \text{elementary charge})$$

i.e. the energy increase $h\Delta\nu$ is just compensated for by the energy loss $e\Delta U_0$.

If negative anode potential $U_0(\nu)$ is plotted against ν , we obtain a straight line with a gradient given by (see (2)):

$$\frac{\Delta U_0}{\Delta\nu} = \frac{h}{e}$$

We obtain the following for a known elementary charge e

$$h = e \frac{\Delta U_0}{\Delta\nu} \quad (3)$$


Fig. 2: Compact arrangement for determining h .

- (1) High-pressure mercury lamp
- (2) Slider
- (3) Converging lens
- (4) Slit
- (5) Imaging lens
- (6) Direct-vision prism
- (7) Mirror
- (8) Measuring line connections
- (9) Threaded guide of the swivelling arm of the photocell
- (10) Converging lens with slit diaphragm
- (11) Photocell
- (8.1) Connection for measuring amplifier
- (8.2) 4-mm socket
- (12) Window and dimming slider

Setting up:

Optical adjustment (Fig. 2):


Preparation: cover the inner side of the Plexi-glass plate of the emission window (12) with white paper. Remove the direct-vision prism (6) and the lens with slit diaphragm (10) from the beam path. Make the necessary electrical connections (Fig. 3). Secure the Hg lamp (1) at a distance of 5 mm in front of the compact arrangement housing (light must be prevented from entering by means of slider (2)) and switch on.

- Form an image of the Hg lamp (1) on the slit (4) with the lens (3) (slit in the center of the holder, adjust the lamp and lens).
- From an image of the slit on the emission window with lens (5) by adjusting lens (5) and, if necessary, mirror (7).
- Insert the direct-vision prism (6) so that the Hg line is on the inside at the emission window with the violet line on the outside (ensure that the black rubber collar is fitted in such a way that no light can pass through).
- Screw the photocell into the socket; place the lens with mounted slit diaphragm in the beam path so that the image of a spectral line is formed on the photo cathode.
- If necessary, repeat point 2.

Apparatus:

| | |
|---|-------------|
| 1 compact arrangement for determining Planck's constant | 558 79 |
| 1 high-pressure mercury lamp | 451 15 |
| 1 universal choke | 451 30 |
| 1 photocell for determining h | 558 77 |
| 1 regulated power supply unit, ≥ 5 V DC, 100 mA | e.g. 522 31 |
| 1 measuring resistor 1 k Ω | 536 14 |
| 1 rheostat 1 k Ω | 537 22 |
| 1 I-measuring amplifier D | 532 00 |
| 1 moving coil instrument P | 531 80 |
| 1 multimeter 1 H-LH | 531 50 |
| 2 connecting leads, blue, 50 cm | 501 26 |
| 4 connecting leads, red, 50 cm | 501 25 |
| 2 connecting leads, blue, 1 m | 501 31 |
| 1 connecting lead, red, 1 m | 501 30 |
| 1 connecting lead, black, 1 m | 501 33 |
| 1 connecting lead, blue, 25 cm | 501 21 |

Fig. 3: Electrical connections

Carrying out the experiment:

Partially darken the room. The slit diaphragm of the lens (10) has a marking pin whose shadow is visible on the emission window.

Determine the negative anode potential U_0 for yellow light.

- To do this, adjust the swivelling arm using the threaded guide (9) so that the shadow of the marking pin is on the yellow spectral line.
- Cover the emission window (12) with the dimming slider if you are not working in a completely darkened room.
- Set U' to 0 by means of the rheostat, and measure the photocurrent I_0 ; then, increase U' until $I' = 0$ (range 10^{-11} A). Read off U_0 .
- Repeat points 1.1-1.3 for the green, turquoise, blue and violet lines. Change the spectral colour by moving the swivelling arm.

Measurement example:

Frequencies of the mercury lines in accordance with Table 1.

Table 1

| Color | yellow | green | turquoise | blue | violet |
|-----------|--------|-------|-----------|------|--------|
| Frequency | 5.19 | 5.49 | 6.88 | 6.08 | 7.41 |

$10^{14} \cdot s^{-1}$


Fig. 4: Frequency ν as a function of negative anode potential U_0 at which the photocurrent I_p becomes 0.

Evaluation and results:

U_0 is plotted against ν in Fig. 4. The frequencies of the Hg spectral lines used are taken from Table 1.

I_0 is plotted against ν as a rough measure of light intensity in Fig. 5. There is clearly no systematic relationship. A comparison with Fig. 4 shows that the kinetic energy of the photoelectrons does not depend on the light intensity, but only on its frequency (the curve in Fig. 4 rises linearly, i.e. the energy of the photoelectrons rises linearly although the intensity does not increase with frequency).

From Fig. 4 we obtain a straight-line gradient

$$\frac{\Delta U_0'}{\Delta \nu} = 4,28 \cdot 10^{-15} \text{ Vs.}$$

With $e = 1.6 \cdot 10^{-19}$ A it follows

$$h = 6,86 \cdot 10^{-34} \text{ Ws}^2.$$

The deviation from the theoretical value, $h =$

$6.625 \cdot 10^{-34} \text{ Ws}^2$, is 3.5 % for this measurement and this photocell. Large deviations are possible due to scattering in the photocell.


Fig. 5: Photocurrent I_p , as a measure of light intensity, as a function of frequency. I_p in arbitrary units.

Notes:

1. The use of potassium as the cathode material is explained by the particularly low work function (alkaline metal with weakly-linked photoelectron).

2. Sources of error:

a) Potassium coating on the anode. The result of this is that this acts as a cathode when the anode ring is illuminated and a photocurrent flows in the "wrong" direction. In particular, a strong "opposing current" occurs when U' is increased through U_0 . This can be remedied by carefully heating out the anode. To do this, connect both anode ring connections to a variable transformer (e.g. 522 20). In a darkened room observe the anode ring and carefully increase the heating voltage. Switch off the heating voltage as soon as the anode ring begins to glow. The whole heating process should not last longer than 15 seconds.

b) The photocell is not clean. The result is a leakage current between anode and cathode which is also measured as an opposing current. This can be remedied by cleaning the photocell with alcohol.

c) The optical surfaces are not clean (dust, fingerprints). The result of this is stray light leading to an excessively high U_0 at low frequencies. This can be remedied by cleaning the optical surfaces, in particular the prism.

The measurement error is mostly less than 10 % if the experiment is carried out carefully.

3. It is principally impossible to measure the work function A_C of the photocathode, in spite of

the view represented in many text books. This is shown by the following observation:

Electrons in the cathode are in a potential well with a depth of A_C (Fig. 6) the latter being measured from the fermi level ϵ_F . In the same way,

the electrons of the anode are also in a potential well with a depth of A_a . If a conducting connection is created between the cathode and anode, the fermi levels adjust so that they are at the same height (Fig. 6).

If a voltage U is applied between anode and cathode, the fermi levels are displaced by $e \cdot U$ with respect to each other. Electrons from the cathode must now overcome the highest point of the potential "mountain" which has been created in order to reach the anode; however, this is given by $A_a + e \cdot U'$.

The anode current becomes zero when

$$h\nu = eU_0 + A_a$$

(Fig. 7).

The often-quoted equation

$$h\nu = eU_0 + A_C$$

assumes that the potential difference $e \cdot U$ is applied between the edges of the potential wells instead of (as is correct) between the fermi levels.

It is not possible to measure A_a without going to a lot of trouble since A_a greatly depends on the crystal orientation and surface accumulation of residual gas atoms.


Fig. 7: Energy relationships when a negative anode potential is applied.


Fig. 6:

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En la figura 5 está graficado I_0 , como medida grosera para la intensidad de luz, contra ν . Apparently no systematic relation exists. A comparison with figure 4 shows that the kinetic energy of the photoelectrons does not depend on the intensity of the light, but only on its frequency (the curve in figure 4 increases linearly, that is, the energy of the photoelectrons increases linearly, although the intensity does not increase with the frequency).

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