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
**Photocell for  $h$ -Determination**  
**and**  
**Housing for the Photocell**

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**1. Description**

The photocell for  $h$ -determination is used to carry out experiments on the photoelectric effect, and in particular — by way of its spectral dependency — to determine Planck's quantum of action.

The **photocell (558 77)** is a potassium-vacuum photocell of large area with an electrically heated anode made of platinum wire. The heater voltage is supplied through an E 14 screw cap, while the photosensitive alkali layer (cathode) is connected up on the opposite side through a metal cap 9 mm in diameter.

**Notes**

- The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
- The data DC ... refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSICS LEAFLETS".  
 Published to date:

	Cat. No.
1st series (72 sheets) . . . . .	599 21
2nd series (52 sheets) . . . . .	599 22
3rd series (49 sheets) . . . . .	599 23
4th series (50 sheets) (translation prepared) . . . . .	599 24

- The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the latest knowledge acquired in all scientific and technical fields.

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The photoelectric saturation current amounts to several  $10^{-7}$  A/lm, and is correspondingly less at low driving potentials.

The **housing for the photocell (558 78)** consists of two parts. The bottom part is mounted on a shaft of 10 mm diameter, and has an E 14 socket to fit the photocell. The position of the screwed-in photocell can be adjusted with the aid of 3 knurled-head screws located at the bottom of the base plate. A two-wire cable provides the anode connections. For connecting up an earth wire, a transverse drill-hole is provided at the top, a centric drill-hole at the bottom of the shaft, both of 4 mm  $\phi$ .

The top part, 58 mm in diameter, can be fitted onto the bottom part, while from the inside a contact is cushioned on the cathode cap of the cell. This connection is highly insulated and leads to a 4 mm socket on the top part of the housing, which on applying slight pressure from above can be turned and taken out. Laterally there is an opening 30 mm in diameter for the admission of light.

## 2. Operation

The photocell is screwed into the E 14 socket on the bottom part of the housing. Subsequently the top part is mounted with its opening in such a position that the cathode can be exposed to light. The apparatus is fitted to a stand, and the 10 mm shaft and one pole of the anode connections are joined to an earth wire. The insulated cathode connection is connected via a sensitive ammeter to a small, controllable D. C. voltage (0... + 2 V against earth.) The demonstration multirange meter (531 86) or the projection type twin moving coil instrument (442 73/74) possesses in combination with the measuring amplifier (532 01) a sufficient sensitivity.

The platinum loop (anode) is heated to vaporize any potassium that has condensed on it from the photosensitive layer. Otherwise this potassium on the anode might cause a reverse electron current if the platinum loop were illuminated, which would falsify the results. But as heating the platinum loop might cause an undue rise of temperature in the photosensitive alkali layer of the photocathode, therefore rendering it useless, the following procedure should be adhered to:

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A source of D. C. or A. C. yielding 2 V should be connected to the cables of the case leading to the E 14 Edison socket, interposing a switch. When switching on, the heating current may be 1 ... 1.2 A, but the switch should only be switched on when the whole arrangement has been set up and the measuring amplifier used to measure the current has been connected up and taken into operation. The platinum loop should be at a potential which is 1 to 2 V below that of the photosensitive layer (reverse potential). The switch for the heating current should be switched on briefly prior to (not during) a determination and then switched off again at once when the moving coil instrument connected to the measuring amplifier begins to indicate a reverse current. This will generally be the case after 1 ... 2 seconds and is the best criterion for a just-sufficient heating of the platinum loop. In no case must the platinum loop be heated over prolonged periods and with somewhat smaller currents, as this would cause the cell to warm up and change the properties of the photosensitive layer.

Heating the platinum loop is required only if—despite cleaning of the housing and cell, and a good optical adjustment—there is a comparatively large counter-current with increasing counter-voltage.

### 3. Hints

The light of the mercury vapour lamp (451 15) filtered with the light filter for the mercury line  $\lambda = 436$  nm (468 32) or the ultraviolet filter  $\lambda = 366$  nm (469 79) respectively, suffices already for an approximate determination of Planck's quantum of action. The school mirror galvanometer (532 10) is used as the indicating instrument.

An exact determination is only possible with spectrally pure light which may be obtained with the mercury vapour lamp (451 15), the direct vision prism (466 05), the adjustable slit (460 14) and appropriate image lenses.

By means of the measuring amplifier (532 01) the moving coil instrument (531 86 or 442 73/74) becomes so sensitive that it is possible to take measurements with wavelengths in the visible range, i. e. 578 nm (yellow), 546 nm (green), 493 nm (green-blue), 436 nm (blue) and 405 nm (violet).

Unsatisfactory measurements may be due to poor insulation of the cell or to a photoelectric effect at the anode. The insulated connection to the housing as well as the photocell must, therefore, be clean at all times. If necessary, they should be cleaned with

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Green's lens tissue (305 00), or with alcohol and subsequent treatment with warm air.

To keep the photoelectric effect at the anode at a minimum, avoid direct exposure of the anode to light. It is recommended to fix a circular diaphragm, 15 mm in diameter, made of black paper, onto the glass bulb of the photocell. Exposure with scattered light from the cathode is unavoidable. With pure platinum, only the mercury resonance line  $\lambda = 254$  nm can release photo-electrons. Occasional heating of the anode therefore improves the result.

#### 4. Experiments

For detailed description on this experiment refer to Leybold Physics leaflets DC 535.215.2; b, 3rd series, sheet 4 Planck's quantum of action  $h$ .

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